



Powerheart AED G3 Plus

The Powerheart® AED G3 Plus is the flagship Cardiac Science automated external defibrillator (AED) that features Rescue Ready®, RescueCoach™ and STAR® biphasic technology to make your rescue as simple as possible. Available in automatic and semi-automatic models, this AED is ideal for lay rescuers to use in public access defibrillation areas such as schools, athletic fields, and workplace environments.

Operations	9390A (fully automatic version) and 9390E (semi-automatic version)
Waveform:	STAR® biphasic truncated exponential
Energy (J) range available:	Escalating variable energy (VE) 95J to 351J
Protocols:	5 energy protocols available
Factory default (nominal):	200VE, 300VE, 300VE
Voice prompts:	RescueCoach voice instructions guide user confidently through rescue process
Text screen:	Displays rescue prompts to guide user through rescue process as well as additional critical rescue information for EMS responders
Visible indicators:	Rescue Ready status indicator, SmartGauge battery status indicator, service indicator, PAD indicator, text display
Audible alerts:	Voice prompt, system alert
Synchronized shock:	Built-in automatic synchronization feature
Pacemaker pulse detection:	Yes
Programmable:	Yes, via MDLink®
Pediatric capability:	Yes
Warranty:	7 years
CPR cadence:	Metronome for compression frequency



Powerheart AED G3 Pro Plus

The Powerheart® AED G3 Pro is a fully equipped automated external defibrillator (AED) that features Rescue Ready® and STAR® biphasic technology and is designed for medical professionals. Featuring a color display, 3-lead ECG monitoring capability and manual defibrillation override, this AED provides you with the flexibility and features you need when saving a life.

Operations	Semi-automatic with manual override
Waveform:	STAR® biphasic truncated exponential
Energy (J) range available:	Escalating variable energy (VE) 95J to 351J
Protocols:	5 energy protocols available
Factory default (nominal):	200VE, 300VE, 300VE
Voice prompts:	Clear, concise voice prompts guide user through the rescue
Visible indicators:	Rescue Ready status indicator
Audible alerts:	Voice prompt, system alert
Synchronized shock:	Built-in automatic synchronization feature
Pacemaker pulse detection:	Yes
Programmable:	Yes, via MDLink®
Pediatric capability:	Yes
Warranty:	7 years
Control buttons:	Shock button and manual override
Display content:	Displays written instructions to guide user through rescue process, SmartGauge battery status indicator, service indicator, pad indicator, text display, ECG display
Display specifications:	3.5 in (8.9 cm) diagonal transreflective TFT display with 320 x 240 pixels (quarter VGA). Resolution is 113.5 dots/in (4.47 dots/mm)

AED Accessories



This **Ready Kit** clips onto the carry bag of the Powerheart AED and provides key tools for responding to sudden cardiac arrest. The kit includes: nitrile gloves, razor, scissors, towel, 4" gauze, one-way filter mask, and carabineer attachment.



This **Intellisense Lithium Battery** is for use in Powerheart AED G3 9300E, 9300A, 9390E, and 9390A devices shipped after April 12, 2004. This non-rechargeable battery comes with a four-year warranty.



The **Adult AED Defibrillation Pads** are non-polarized so that either pad can be placed in either location, simplifying the rescue. These defibrillation pads can be used with all Powerheart Automated External Defibrillators.



Pediatric Defibrillation Pads
These reduced-energy child/pediatric defibrillation pads work with all 9200 and 9300 series AEDs, enabling rescuers to provide defibrillation therapy to children who are less than 8 years of age or who weigh less than 55 pounds



Carry Bag

This semi-rigid carry bag provides additional protection for your Powerheart G3 AED and makes it easy to grab the unit in a hurry. The large back pouch is perfect for storing extra electrodes and the shoulder strap is easily adjustable. In addition, the Rescue Ready indicator and pads' expiration date are visible at all times. Fits all 9300 series AEDs.



Wall Storage Sleeve

This AED wall mounted storage sleeve is an economical and attractive AED holder which keeps your Automated External Defibrillator in an easily accessible location.



Wall Cabinet

This well constructed surface mounted AED wall cabinet comes with an alarm and highly visible strobe light that are triggered when the door is opened — making it widely-known that an emergency is in progress. This alarm can be wired into a building's security system and can be enabled or disabled with a key that is provided. The AED sits inside the white cabinet with a see through door.



Recue Backpack

This AED Rescue Backpack is a great storage and carry solution for users who take their Powerheart AED with them on the go, including bicycle medic teams, search and rescue teams, and scout troops. The backpack stores the AED in an easily accessible outer pouch. The spacious main compartment is filled with pouches ideally suited to organize, store, and carry other essential rescue supplies.



Waterproof Carrying Case

This hard-sided waterproof carry case is a great storage solution for rugged environments. The Powerheart AED fits safely and snugly inside while the durable outer shell provides additional protection against bumps, drops and the elements.